

The Surface Morphology of Normal Human Leukocytes by Chilled Scanning Electron Microscopy

Akiko Adachi, Shigeru Sato, Yoshihiro Sasaki and Zenya Naito
Division of Morphological and Biomolecular Research, Nippon Medical School




Fig. 1




Fig. 2

In conventional scanning electron microscopy (SEM), granulocytes are difficult to distinguish from monocytes and lymphocytes. Therefore, we tried a new method of chilled SEM¹. Chilled SEM, a low-vacuum SEM with a cooling stage, provides images through back-scatter electrons. In comparison with specimens for conventional SEM, specimens for chilled SEM are not completely dehydrated in ethanol and require no critical-point drying and no metal coating¹. Chilled SEM images provide information beneath the surfaces of uncoated specimens^{2,3}. Thus, chilled SEM appears able to distinguish different leukocytes. Granules are observed in neutrophils, eosinophils, and basophils. Lymphocytes have no granules. Monocytes have well-developed ridge-like ruffles that other leukocytes do not.

Correspondence to Akiko Adachi, Division of Morphological and Biomolecular Research, Nippon Medical School, 1-1-5 Sendagi, Bunkyo-ku, Tokyo 113-8602, Japan




Fig. 3




Fig. 4

- Fig. 1** Comparison of conventional SEM and chilled SEM for lymphocytes. **a:** The surface of cells have many microvilli on conventional SEM. **b:** Many smaller ridge-like ruffles are seen on chilled SEM. Circulating disk-shaped red blood cells (RBCs) have the same surface morphology with both conventional SEM and chilled SEM.
- Fig. 2** Ultrastructure of monocytes. **a:** A well-developed cytoplasmic projection (**arrows**) is seen on transmission electron microscopy. Pl, platelet. **b:** Well-developed ridge-like ruffles (**arrows**) are seen with chilled SEM.
- Fig. 3** Ultrastructure of neutrophils. **a:** Cytoplasmic projection and many small granules (**arrows**) are seen with transmission electron microscopy. **b:** Smaller ridge-like ruffles and a few small granules (**long arrows**) in neutrophils are seen with chilled SEM. A few small granules (**short arrows**) are also seen in platelets (PI).
- Fig. 4** Ultrastructure of eosinophils. **a:** Cytoplasmic projection and many big granules (**arrows**) are seen with transmission electron microscopy. **b:** Big granules (**arrows**) and smaller ridge-like ruffles are seen with chilled SEM.

References

1. Sasaki Y, Sato S, Adachi A, Dan Y, Nishimura M: Use of oolong tea extract staining soft-tissue specimens in low-vacuum scanning electron microscope with a cooling stage. *Med Electron Microsc* 2001; 34: 254-257.
2. De Harven E: Scanning electron microscopy in the backscattered electron imaging (BEI) mode: Application to clinical hematology. *Ultrastruct Pathol* 1987; 11: 711-721.
3. Sato S, Matsui H, Sasaki Y, et al.: The efficiency of X-ray microanalysis in low-vacuum scanning electron microscope: deposition of calcium on the surface of implanted hydrogel intraocular lens (IOL). *J Submicrosc Cytol Pathol* 2006; 38: 1-4.