-Photogravure-

Chondrogenic and Osteogenic Differentiation of Adipose-derived Stem Cells Isolated from GFP Transgenic Mice

Rei Ogawa¹, Hiroshi Mizuno¹, Hiko Hyakusoku¹, Atsushi Watanabe², Makoto Migita² and Takashi Shimada² ¹Department of Plastic and Reconstructive Surgery, Nippon Medical School ²Department of Biochemistry and Molecular Biology, Nippon Medical School



A. Light microscopic images

B. Fluorescence microscopic images

Fig. 1



A. Light microscopic images



B. Fluorescence microscopic images



C. Alcian blue staining

Correspondence to Rei Ogawa, MD, Department of Plastic and Reconstructive Surgery, Nippon Medical School, 1–1–5 Sendagi, Bunkyo-ku, Tokyo 113–8603. Japan. r.ogawa@nms.ac.jp Journal Website (http://www.nms.ac.jp/jnms/)

J Nippon Med Sch 2004; 71(4)



A. Light microscopic images

C. Alkaline phosphatase staining



B. Fluorescence microscopic images



D. von Kossa staining

Fig. 3

Abstract

Recent studies suggest that human adipose tissue contain pluripotent cells similar to bone marrow-derived stromal cells (BSCs)¹⁻³. Taking advantage of homogeneously marked cells from green fluorescent protein (GFP) transgenic mice⁴, we have previously demonstrated that BSCs differentiate into a variety of cell lineages both *in vitro* and *in vivo*⁵. In the present study, we extend this approach to characterize adipose-Integes both *in vitro* and *in vivo*². In the present study, we extend this approach to characterize aupose-derived stromal cells (ASCs)⁶. These cells derived from human are sometimes called processed lipoaspirate (PLA) cells. ASCs were prepared from inguinal fat pads of GFP transgenic mice after extensive washing with PBS and treatment with collagenase. After the primary culture in control medium (DMEM + 10% FBS), the cells were incubated in either chondrogenic medium (DMEM + 1% FBS + insulin + ascorbate 2-phosphate + TGF-beta 1) or osteogenic medium (DMEM + 10% FBS + dexamethasone + ascorbate-2-phosphate + beta-glycerophosphate) for two to four weeks. Chondrogenic differentiation was assessed by Alcian blue staining. while osteogenic differentiation was by von Kossa and Alkaline phosphatase staining. ASCs incubated in chondrogenic medium induced Alcian blue positive cells. Incubation with osteogenic medium became positive for von Kossa and Alkaline phosphatase staining. No osteochondrogenic differentiation was observed in cells incubated with control medium. This cell population can be easily identified through fluorescence microscope, it should be an ideal source of ASCs for further experiments of stem cell biology and tissue engineering.

References

- Zuk PA, Zhu M, Mizuno H, et al: Multilineage cells from human adipose tissue: Implications for cellbased
- therapies. Tissue Eng 2001; 7: 211–228. Mizuno H, Zuk PA, Zhu M, et al: Myogenic differentiation by human processed lipoaspirate cells. Plast Reconstr Surg 2002; 109: 199.
- Zuk PA, Zhu M, Ashijan P, et al: Human adipose tissue is a source of multipotent stem cells. Mol Cell Biol 2002; 13: 4279–4295. 3.
- Okabe M, Ikawa M, Kominami K, et al: 'Green mice' as a source of ubiquitous green cells. FEBS Lett 4. 1997; 407: 313-319.
- Hayakawa J, Migita M, Ueda T, et al: Generation of a chimeric mouse reconstituted with green fluorescent protein-positive bone marrow cells: a useful model for studying the behavior of bone marrow cells in regeneration in vivo. Int J Hematol 2003; 77: 456–462. Ogawa R, Mizuno H, Watanabe A, et al: Osteogenic and Chondrogenic Differentiation by Adipose-derived Stem Cells Harvested from GFP Transgenic Mice. BBRC 2004; 313: 866–872. 5.

Fig. 1 Primary culture of ASCs-cultured for 4 days in control media-

The GFP transgenic mouse ASCs appeared as a fibroblast-like cells (A). GFP positive cells can be easily identified through fluorescence microscopy. **Fig. 2** Two weeks of "culture" in chondrogenic media GFP transgenic mouse ASCs in the micromass culture supplemented with chondrogenic media formed

small spheroids (A, B). The nodules stained positively for Alcian blue staining (C).

Fig. 3 Four weeks of culture in osteogenic media GFP transgenic mouse ASCs in osteogenic media changed fibroblastic appearance to a rounder, more cuboidal shape (A, B). These cells stained positively in Alkaline phosphatase staining (C), and extracellular matrix stained positively in von Kossa staining (D).